

En Novative Technologies, Inc.

En Core®



### 5 Gram Disposable En Core® Sampler



Delivers approximately 5 grams of soil. Use in accordance with SW 846 Method 5035. Consult your regulatory agency for detection limit requirements to help determine which size En Core Sampler to use. Works with Standard T-Handle.

### 25 Gram Disposable En Core® Sampler



Delivers approximately 25 grams of soil. Use in accordance with SW 846 Method 5035. Consult your regulatory agency for detection limit requirements to help determine which size En Core Sampler to use. Works with Standard T-Handle.

### Standard En Core® T-Handle



Attach either a 5 gram or 25 gram En Core disposable sampler to this Standard T-Handle and push the sampler into the soil for a measured soil sample. T-Handle measures ~4.5" long. T-Handle extensions are available. Call for details.

### En Core® Extrusion Tool

When either a 5 gram or 25 gram disposable En Core Sampler arrives at the laboratory, attach it to the Extrusion Tool, engage the plunger, and push the sample into an appropriate vial.



*Special Offer!*

## Disposable En Core Sampler Starter Kit

Kit contains 1 T-Handle, 12 En Core Samplers, 1 Demo Sampler, and a Certificate good for \$100 toward the purchase of 3 cases of En Core Samplers. Please call Toll Free 1-888-411-0757 for details and pricing.



Check us out online  
at <http://www.ennovativetech.com>

En Novative Technologies, Inc.  
1241 Bellevue Street • Green Bay, WI 54302  
Email: [info@ennovativetech.com](mailto:info@ennovativetech.com)  
Phone: 920-465-3960 • 920-465-3963  
Toll Free: 888-411-0757